



EyeSpy 20/20™

Quick Start Guide

INSTALL



1 GO TO VisionQuest2020.org, click on EyeSpy 20/20 menu, & select "Install". Follow on-screen instructions.

2 CLICK RUN to begin the installation process.

3 LAUNCH EYESPY 20/20 by double-clicking the EyeSpy 20/20 icon on the desktop and click OK when asked.

4 IMPORTANT - Slide the lever all the way right to the 'Unlimited' setting and click CLOSE.



5 ENTER your 16-digit activation code and click the "Submit" button.

6 ACCEPT the User Agreement. On the next screen, click OK to quit the EyeSpy 20/20 program.

7 RESTART EYESPY 20/20 by double-clicking the EyeSpy 20/20 icon on the desktop and click OK when asked.

8 CLICK continue when your school appears on the screen.

9 CALIBRATE when prompted by aligning a \$1 bill on the screen and clicking the lower right corner of the bill. When the box is the correct size, click SAVE.

Welcome to EyeSpy 20/20

SET-UP

1

PLACE COMPUTERS on a table and plug in to electrical outlet.



2

ADJUST SCREEN of the laptop computers to a 90 degree angle to ensure there is no glare.



3

POSITION additional table/desk with a chair 10 feet from the laptops. Use the 10 foot mouse extension cable to measure distance between screen and the child.



4

CONNECT MOUSE

Plug one end of the mouse extension cable into the computer and plug the other end into the mouse on the table where the children will be seated.

5

LEAVE 2 FEET OF SPACE

between each laptop and children's seats allowing screening proctors to move about easily without blocking the view of the children playing the game.



6

EQUIP each station with a mouse, stereo glasses, and occluding glasses.

7

TURN ON COMPUTERS

and start the EyeSpy 20/20 software. An internet connection while screening is NOT required.



SCREEN

1 ENTER STUDENT INFO

before conducting screenings. Import data from preloaded student rosters or manually enter data.

EyeSpy 20/20™ Click here for demo mode

VisionQuest 20/20 Demonstration School

Proctor Name Student ID

Student FIRST Name Student LAST Name

DOB (mm/dd/yyyy) Test Color Vision

No answer provided

Ethnicity (select from list)

Grade (select from list) Teacher Name / HR

Female Male English Spanish

Gender Default Report

None Wearing Forgot Lost / Broken

Prescription glasses or contacts

SELECT student from list How to vision screen?

View License * Require Student ID

Game Volume Check Volume

Test Distance 10 Calibration

Reports School Password *

Input Near VA No animation Acquire Exact VA

Vision Reports Detailed Flies/Fall Overst only

View License * Require Student ID

Game Volume Check Volume

Test Distance 10 Calibration

Password Verify password LOCK

BEGIN screening

2:24:05 PM - Thursday, January 5, 2017 - Ver 7.35 Y001K258A19A - U.S. Patent No. 8,808,267

2 ENTER PROCTOR NAME

in the upper left box on the screen. This only needs to be entered once per screening session.

3 VERIFY status of child's glasses.

None Wearing Forgot Lost / Broken

Prescription glasses or contacts

4 BEGIN SCREENING

by clicking the "Begin Screening" button. Screening begins when child chooses a character using the mouse & clicking.

5 CHILD PLAYS GAME



6 FOLLOW on-screen prompts.

ATTENTION PROCTOR !

Make sure student is not wearing an eye patch or colored glasses at this time.

Press Enter when ready.

7 REPORT

appears at end of game. Data automatically saved. Press "S" key to see report in Spanish. Press 'Enter' key to print the report.

IMPORTANT - DO NOT DISCARD !

VISION SCREENING REPORT PREPARED EXCLUSIVELY FOR

Age 8 - DOB 07/28/2008

Teacher: Orta C - Grade: 0

Printed: 01:13:28 PM - Thursday, November 3, 2016 by

VisionQuest 20/20 (www.visionquest2020.org)

in partnership with

OVERALL VISION EVALUATION (wearing glasses)

This vision screening suggests your child may have a vision problem.

If your child is not currently under the care of an eye doctor, a professional eye examination is recommended.

Eye Professional Use Only

Does the child have normal vision without glasses? Y / N

Does the child have normal vision with glasses? Y / N

Are there any concerns about vision? Y / N

Does the child have normal vision with glasses, which would enable safe school performance? Y / N

Does the child have normal vision with glasses, which would enable safe school performance, at school? Y / N

Follow-up appointment: _____

Printed Student's Name: _____

Eye S. One above the S

ATTENTION! ARIZONA TEST PAPERS

YOU MAY DISCARD THIS REPORT AFTER YOU HAVE PRINTED IT. THIS REPORT IS NOT VALID FOR ANY OTHER PURPOSE.

For more information please contact: Barbara Turner | Health Office @

visionquest2020.org | 800-368-3688 | 800-368-3688

Screening is complete !

Press the Enter key to print the screening report

OR

Press the spacebar to screen the next student

Press the S key to show the report in Spanish

REPORTS

1

CLICK REPORTS from the main screen of EyeSpy 20/20.



4

SORT & SEARCH

The list of reports can be sorted by pressing the SORT button above the column by which sorting is desired.

Search reports by entering criteria in one or more search fields above the column(s) by which search is desired. The resulting list of reports will contain only records matching the search criteria.

Sorts and Searches can be combined to display reports based on specific criteria.

2

ENTER DATE RANGE of the reports you wish to see or click the starting month to see screenings through today. Click the appropriate SHOW button and the reports will appear.



5

PRINT & EXPORT

the displayed vision screening records via the buttons at the bottom of the screen. Summary and individual reports are accessed by the appropriately named buttons. Vision screening data can also be exported to a Microsoft Excel file.

3

SELECT CHILD from the list of reports. Double-click the child's name to display the child's vision screening report.

